

16. BUFFER PROPERTIES

Initial Question

Buffers are solutions that are resistant to changes in their pH when acids or bases are added. For example, human blood contains the bicarbonate ion. This ion can accept hydrogen ions to remove excess acidity in the blood or can donate hydrogen ions to remove alkalinity in the blood. Once the bicarbonate ions are used up, blood can rapidly become either too acidic or too basic. In other words, the bicarbonate buffer system in blood has a limited capacity.

How are buffers made, and what determines their capacity?

Materials and Equipment

Model 1

- Data collection system
- pH sensor
- Analytical balance
- Volumetric flask, 100-mL or 250-mL¹
- Beakers (2), glass, 50-mL
- Sodium acetate (NaCH₃COO), about 1.0 g²
- Ammonium chloride (NH₄Cl), about 1.0 g²
- 0.3 M Acetic acid (CH₃COOH), 100 mL²
- 0.3 M Ammonia (NH₃), 100 mL²
- Distilled water, 150 mL²
- Marking pen

¹The volume of the flask depends on the buffer assigned from the Model 1 Data Table.

²The volume and mass needed depends on the buffer assigned from the Model 1 Data Table.

Model 2

- Universal indicator, 3 drops
- Beral pipets (2)
- Test tubes (3), 20 mm × 150 mm, glass, 25-mL
- 0.10 M Hydrochloric acid (HCl), 20 mL
- 0.10 M Sodium hydroxide (NaOH), 20 mL
- Buffer solution from Model 1
- Distilled and deionized water

Applying Your Knowledge

- Data collection system
- pH sensor
- Beakers (2), 50-ml
- Volumetric flask, 100-ml
- Two of the following, to create 100 mL of buffer:
 - 0.3 M Acetic acid (CH₃COOH)
 - 0.3 M Sodium acetate (NaCH₃COO)
 - 0.3 M Sodium phosphate dibasic (Na₂HPO₄)
 - 0.3 M Sodium phosphate monobasic (NaH₂PO₄)
 - 0.3 M Ammonia (NH₃)
 - 0.3 M Ammonium chloride (NH₄Cl)
 - 0.3 M Potassium phosphate (K₃PO₄)
 - 0.3 M Phosphoric acid (H₃PO₄)
- 6 M Sodium hydroxide (NaOH), 5 drops¹
- 6 M Hydrochloric acid (HCl), 5 drops¹
- Stirring rod

¹Your teacher will add the 5 drops of NaOH and HCl to test your buffer.

Safety

Add these important safety precautions to your normal laboratory procedures:

- Wear your goggles.
- This lab uses strong acids and bases. If you come in contact with either of them, flush the area with plenty of water.

Getting Your Brain in Gear

1. For each solution described below, determine if it acts as a buffer. Explain your reasoning.

a. 0.50 mole of HCl and 0.50 mole of NaCl in 250 mL of water

b. 0.50 mole of HNO_2 and 0.50 mole of NaNO_2 in 250 mL of water

c. 1.00 mole of CH_3COOH and 0.50 mole NaOH in 250 mL of water

2. Consider a buffer made of an equal number of moles of HF and NaF.

a. Write a net ionic equation to illustrate how the buffer solution reacts to keep the pH stable when strong acid is added.

b. Write a net ionic equation to illustrate how the buffer solution reacts to keep the pH stable when strong base is added.

MODEL 1

Building Model 1 – Buffer Preparation

1. After your instructor assigns you a buffer from Model 1 to prepare, refer to the Model 1 Data Table to determine the chemicals and the molarities you will need for your buffer.

Solution assigned: _____

2. What component(s) of your assigned buffer system are available (refer to the Materials and Equipment list for Model 1)?

3. Calculate the number of moles of weak acid and conjugate base, or weak base and conjugate acid, you will need in the final solution to make the buffer you have been assigned.

4. Calculate the mass of the solid and volume of solution you will need to make your assigned buffer.

5. Obtain the appropriate size volumetric flask.

6. Fill the volumetric flask halfway with distilled water.

7. Obtain _____ grams of solid _____ and add it to the volumetric flask.

NOTE: Do not fill the flask with water.

8. Use a graduated cylinder to measure _____ milliliters of _____ and add it to the volumetric flask. Swirl the flask to mix the solution.

9. Fill the volumetric flask to the mark with distilled water. Mix the solution thoroughly by inverting the flask several times.

10. Rinse a 50-mL beaker three times with distilled water.

11. Rinse the 50-mL beaker twice with small portions of your buffer solution.

12. Pour about 25 mL of your buffer solution into the 50-mL beaker and measure the pH of your solution. Record the pH in the Model 1 Data Table.

13. Why is it necessary to rinse the 50-mL beaker with distilled water, and then with the solution you want to test?
-
-
-

14. Label and save your remaining buffer for building Model 2. Share your pH data with your classmates to complete the Model 1 Data Table.

Model 1 – Buffer Preparation

Table 1: Model 1 Data Table—pH of buffers of different volumes and concentrations

Solution Number	Molarity of Weak Base	pKa	Molarity of Conjugate Acid	Ratio of Weak Base to Conjugate [OH] / [B ⁺]	Total Volume of Buffer Solution	pH of Buffer Solution	pK _a - pH
1	0.05 M NH ₃		0.05 M NH ₄ ⁺		100 mL		
2	0.10 M NH ₃		0.10 M NH ₄ ⁺		100 mL		
3	0.10 M NH ₃		0.10 M NH ₄ ⁺		250 mL		
4	0.15 M NH ₃		0.15 M NH ₄ ⁺		100 mL		
5	0.15 M NH ₃		0.05 M NH ₄ ⁺		100 mL		
6	0.05 M NH ₃		0.15 M NH ₄ ⁺		100 mL		
Solution Number	Molarity of Weak Acid	pKa	Molarity of Conjugate Base	Ratio of Weak Acid to Conjugate [HA] / [A ⁻]	Total Volume of Buffer Solution	pH of Buffer Solution	pK _a - pH
7	0.05 M CH ₃ COOH		0.05 M CH ₃ COO ⁻		100 mL		
8	0.10 M CH ₃ COOH		0.10 M CH ₃ COO ⁻		100 mL		
9	0.10 M CH ₃ COOH		0.10 M CH ₃ COO ⁻		250 mL		
10	0.15 M CH ₃ COOH		0.15 M CH ₃ COO ⁻		100 mL		
11	0.15 M CH ₃ COOH		0.05 M CH ₃ COO ⁻		100 mL		
12	0.05 M CH ₃ COOH		0.15 M CH ₃ COO ⁻		100 mL		

Analyzing Model 1 – Buffer Preparation

15. Use a reference source to look up the K_a and determine the pK_a of the weak acid or conjugate acid involved in each of the solutions in the Model 1 Data Table.

16. Calculate the weak acid or weak base to conjugate ratio for each buffer solution in Model 1.

- a. Which ratios of weak acid or weak base to the conjugate solution have the pH of the buffer the closest to the pK_a ?

17. Compare the buffers having the same molarity of components, but with different total volumes. Does the volume of the buffer system prepared affect the pH of the resulting solution? Provide specific examples from Model 1 to support your answer.

18. How would you respond to a person who mistakenly said “Buffers are made so that the solution remains neutral”? Provide specific examples from Model 1 to support your answer.

19. Suppose you wanted to make a buffer solution that would keep the pH near 3. Which of the following mixtures would suffice?

Solution A 0.10 M lactic acid and 0.10 M sodium lactate

Solution B 0.10 M hypochlorous acid 0.10 M sodium hypochlorite

Solution C 0.10 M benzoic acid and 0.10 M sodium benzoate

MODEL 2**Building Model 2 – Buffer Capacity**

1. Rinse three test tubes several times with distilled water. Assign one pipet for the deionized water, one for the buffer, and one for the universal indicator.
2. Place 10 drops of deionized water in a test tube. Add one drop of universal indicator to the water and note the initial color. Add 0.10 M HCl solution drop-wise, mixing between drops, until the solution changes to red-orange ($\text{pH} \leq 3$). In Model 2, record the number of drops of HCl solution added to make this color change occur.
3. Place 10 drops of your buffer solution in a test tube. Add one drop of universal indicator to the solution and note the initial color. Add 0.10 M HCl solution drop-wise, mixing between drops, until the solution changes to red-orange ($\text{pH} \leq 3$). Record the number of drops of HCl solution added to make this color change occur.

NOTE: If your solution contains CH_3COOH and CH_3COO^- , then don't carry out this step.

4. Place 10 drops of deionized water in a test tube. Add one drop of universal indicator to the water, and note the initial color. Add 0.10 M NaOH solution drop-wise, mixing between drops, until the solution changes to purple ($\text{pH} \geq 11$). Record the number of drops of NaOH solution added to make this color change occur.
5. Place 10 drops of your buffer solution in a test tube. Add one drop of universal indicator to the solution, and note the initial color. Add 0.10 M NaOH solution drop-wise, mixing between drops, until the solution changes to purple ($\text{pH} \geq 11$). Record the number of drops of NaOH solution added to make this color change occur.

NOTE: If your solution contains NH_4^+ and NH_3 , then don't carry out this step.

6. Share your data with your classmates to complete the Model 2 Data Table.

Model 2 – Buffer Capacity

Drops of HCl needed in distilled water to turn the color red-orange: _____

Drops of NaOH needed in distilled water to turn the color purple: _____

Table 2: Model 2 Data Table—Comparing buffer capacity

Solution Number	Molarity of Weak Base	Molarity of Conjugate Acid	Total Volume of Buffer Solution	Drops of HCl Needed	Drops of NaOH Needed
1	0.05 M NH ₃	0.05 M NH ₄ ⁺	100 mL		
2	0.10 M NH ₃	0.10 M NH ₄ ⁺	100 mL		
3	0.10 M NH ₃	0.10 M NH ₄ ⁺	250 mL		
4	0.15 M NH ₃	0.15 M NH ₄ ⁺	100 mL		
5	0.15 M NH ₃	0.05 M NH ₄ ⁺	100 mL		
6	0.05 M NH ₃	0.15 M NH ₄ ⁺	100 mL		
Solution Number	Molarity of Weak Acid	Molarity of Conjugate Base	Total Volume of Buffer Solution	Drops of HCl Needed	Drops of NaOH Needed
7	0.05 M CH ₃ COOH	0.05 M CH ₃ COO ⁻	100 mL		
8	0.10 M CH ₃ COOH	0.10 M CH ₃ COO ⁻	100 mL		
9	0.10 M CH ₃ COOH	0.10 M CH ₃ COO ⁻	250 mL		
10	0.15 M CH ₃ COOH	0.15 M CH ₃ COO ⁻	100 mL		
11	0.15 M CH ₃ COOH	0.05 M CH ₃ COO ⁻	100 mL		
12	0.05 M CH ₃ COOH	0.15 M CH ₃ COO ⁻	100 mL		

Analyzing Model 2 – Buffer Capacity

The *buffer capacity* of a solution is related to the number of moles of acid and base that the solution can neutralize without a significant change in pH.

For each set of solutions shown below, consider the similarities and differences in their components. Refer to both the Model 1 Data Table and the Model 2 Data Table.

7. Compare the buffer volume, molarity, and pH:

a. What is the molarity ratio of weak base to conjugate acid in solutions 2 and 3?

b. Was there a change in pH between solutions 2 and 3?

c. What is the molarity ratio of weak acid to conjugate base in solutions 8 and 9?

d. Was there a change in pH between solutions 8 and 9?

e. For solutions 2 and 3, and solutions 8 and 9, was there a difference in the number of drops of HCl or NaOH needed to obtain pH 1 or pH 14?

f. Considering the previous questions, what can be concluded about changing the volume of the buffer system without changing the ratio of weak acid to conjugate base, or weak base to conjugate acid, and its effect on pH?

8. The following questions compare the number of moles of each buffer system component and the buffer capacity in solutions of equal volume:

a. What is the ratio of the number of moles of weak base between solutions 1 and 2 and between solutions 1 and 4?

b. Was there a change in pH between solutions 1, 2 and 4?

c. What is the ratio of the number of moles of weak acid between solutions 7 and 8 and between 7 and 10?

d. Was there a change in pH between solutions 7, 8, and 10?

e. For solutions 1, 2, and 4, and solutions 7, 8, and 10, was there a difference in the number of drops of HCl or NaOH needed to obtain pH 1 or pH 14?

f. Considering the previous questions, what can be concluded about buffer capacity and concentration? Give evidence from your data to support your answer.

9. Do buffer systems stop the pH from changing? Give evidence from your data to support your answer.

10. Do buffers with greater concentrations have greater buffer capacity? Give evidence from your data to support your answer.

11. If you wanted to make a buffer of pH 9.5 with a buffer capacity greater than that in Solution 1 of Model 2, what would have to change in the buffer system?

Connecting to Theory

Buffers are made from a weak acid or base and their respective conjugate. Two nutritional labels are shown below. One is from Monster Energy® drink and the other is from Mountain Dew®. By inspecting each label, you should be able to identify an acid and its conjugate. Together these form a buffer solution.

Remember, when acids are made from polyatomic acids, their ending is changed to “ic.” The polyatomic ion of the conjugate will have an “ate” ending.

Nutrition Facts

Monster MW2 Assault Energy Drink

Serving Size - 8 fl. oz.

Servings per container - 2

Amount Per Serving

Calories - 100

Total Carb - 25g

Sugars - 25g

Vitamin B2 - 1.7mg

Vitamin B3 - 20mg

Vitamin B6 - 2mg

Vitamin B12 - 6mcg

Sodium - 1000mg

Taurine - 1000mg

Panax Ginseng - 200mg

Energy Blend - 2500mg

Ingredients: Carbonated Water, Sucrose, Glucose, Citric Acid, Sodium Citrate, Taurine, Natural Flavors, Panax Ginseng Root Extract, Phosphoric Acid, L-Carantine, Caffeine, Sorbic Acid, Benzoic Acid, Caramel Color, Niacinamide, Sucralose, Sodium Chloride, Glucuronolactone, Inositol, Gurana Seed Extract, Certified Color, Pyrodozine Hydrochloride, Riboflavin, Maltodextrin, Cyanocobalamin

Consume Responsibly - Limit 3 cans per day. Not recommended for children, pregnant women, or people sensitive to caffeine.

Nutrition Facts

Halo 3 Mountain Dew Game Fuel

Serving Size - 1 Can

Amount Per Serving

Calories - 170

Total Carb - 0g

Sodium - 65mg

Total Carbohydrates - 46g

Sugars - 46g

Protein - 0g

Caffeine Content: 73mg/12 fl. oz.

Ingredients: Carbonated Water, High Fructose Corn Syrup, Citric Acid, Natural Flavors, Sodium Benzoate, Gum Arabic, Caffeine, Sodium Citrate, Yellow 5, Glycerol Ester of Wood Rosin, Calcium Disodium Edta, Yellow 6, Red 40, Brominated Vegetable Oil

Halo 3 Mountain Dew Game Fuel is not a significant source of other nutrients.

Applying Your Knowledge – Making Your Own Buffer System

Listed below are solutions that can be combined to make buffers. Your teacher will assign you a pH. Your job will be to make 100 mL of a buffer solution with that pH, using these components.

After your teacher measures the initial pH, you will divide your solution between two beakers.

Your teacher will add 5 drops of 6.0 M HCl to one beaker and measure the pH to determine if it has changed. Five drops of 6.0 M NaOH will be added to the other beaker and similar measurements will be made.

Your first task is to find all of the conjugate relationships in the list below. Show each relationship and the information needed to calculate the pH for each. You need to be ready to prepare the solution that matches the assigned pH.

0.3 M Ammonium chloride (NH_4Cl)

0.3 M Potassium phosphate (K_3PO_4)

0.3 M Sodium acetate (NaCH_3COO)

0.3 M Sodium phosphate monobasic (NaH_2PO_4)

0.3 M Sodium phosphate dibasic (Na_2HPO_4)

0.3 M Acetic acid (CH_3COOH)

0.3 M Phosphoric acid (H_3PO_4)

0.3 M Ammonia (NH_3)