




CONNECT



- 1 Download & install *MatchGraph!* from pasco.com/matchgraph
- 2 Launch *MatchGraph!*
- 3 **BLUETOOTH** – Turn your sensor on and connect it
USB – Connect wireless sensor via USB,
or PASPORT sensor via interface (e.g., 550, USB Link)

PASCO tech support 877.373.0300



MATCH



1. Enter student name

Position

2. Choose Position or Velocity



3. Choose profile



4. Record data



Take Journal Snapshot



Show Journal Snapshots



Share / Export Data

pasco.com/matchgraph/

for more *MatchGraph!*[™]
instructions, ideas, and exercises

